

PA-45

Glass Block, Rectangular, Model PA-45, high quality glass with optically worked surfaces. All sides polished and angles accurately worked. Size : 114 x 63 x 19 mm.



PA-45D


Glass block, semi-circular, **Model PA-45D** high quality clear glass. diam. 90mm, thickness 15mm.



Model PA-45R

Model PA-45R Ray Box has three openings, two of which can be used simultaneously by two sets of pupils. The two side openings each have two slots for accommodating blanking plates. A front opening provides seven grooves for accessories and lenses. Supplied with 1m of flexible plastic coated wires with two 4mm banana plugs, a 12V 20W halogen lamp, slit plate with single and triple slit and 2 blanking plates.

Accessories set for PA-45R Ray box



Model PA-45R-1

Colour mixing set with two 50 x 50mm hinged mirrors and one each red, green and blue filters in specially designed off-set positions to enable colour mixing for the ray box.



Model PA-45R-2

Un-mounted acetate sheet, 75 x 50mm, comprising three primaries and three secondaries. Light blue, dark blue primary, yellow, red primary, pale red and dark green primary. Set of 6,



Model PA-45R-3

Comprising one each plano-convex and plano-concave cylindrical lenses, focal lengths 150 and 300mm respectively. Size : 40 x 22mm (W. x H.)