

**Pulleys**


These economical pulleys are perfect for students use. They come with steel frames, brass axles, wide hubs and rigid plastic free-running sheaves. Sheaves on single axle pulleys are 50 mm in diameter. All pulleys feature steel hooks at each end.

- Model PA-80** Single Sheave
- Model PA-81** Double Sheave
- Model PA-80A** Triple Sheave



PA-80A, PA-81, PA-80

PA-81A



**Model PA-81A** Double in line,  
Wheels diameter are 30 & 40 mm.


**Model PA-81B** Triple in line,  
Wheels diameter are 30, 40 & 50 mm.

PA-81B



PA-83-2, PA-83, PA-82

- Model PA-82** Aluminium pulley, single sheave, mounted in metal frame, wheel diameter 50 mm.
- Model PA-83** Aluminium pulley, double sheave, mounted in metal frame, wheel diameter 50 mm.
- Model PA-83-2** Aluminium pulley, triple sheave, mounted in metal frame, wheel diameter 50 mm



PA-83-3

**Model PA-83-3** Aluminium pulley, double in line,  
mounted in metal frame, wheel diameters are 40 & 50 mm.

**Model PA-83-1** Aluminium pulley, triple in line, mounted in metal frame, wheel diameters are 30, 40 & 50 mm.

PA-83-1