

SLOTTED WEIGHT & HANGER

Slotted Weights, economical solid stainless steel, flat, with large central hole and small slot to prevent masses sliding off hanger. For general lab use at the junior and senior school level. Each marked with its weight.

Model

PA-85AS	1 gm
PA-85BS	2 gm
PA-85CS	5 gm
PA-86S	10 gm
PA-87S	20 gm
PA-88S	50 gm
PA-89S	100 gm



PA-85AS to PA-89S



PA-90S to PA-92S

Model

PA-90S	10 gm
PA-91S	20 gm
PA-91AS	50 gm
PA-92S	100 gm

Weight hanger for use with slotted weights, stainless steel. The hanger marked with its weight. These weight hangers will only work with the slotted weights as listed and might not work with any other brand.

HOOK WEIGHT SET

Model PA-92A, for general laboratory use and in particular, with trip scales. Made of solid metal with hook at top and recessed bottom with crossbar. Value is marked on each weight. Set of 9 in a plastic block (1 x 10 gm, 2 x 20 gm, 1 x 50 gm, 1 x 100 gm, 2 x 200 gm, 1 x 500 gm and 1 x 1000 gm). Total 2,100 gram.



PA-92A



PA-92A-1 to PA-92A-7

Single Hook Weights, replacement for model PA-92A

Model

PA-92A-1	10 gm
PA-92A-2	20 gm
PA-92A-3	50 gm
PA-92A-4	100 gm
PA-92A-5	200 gm
PA-92A-6	500 gm
PA-92A-7	1000 gm

METRIC WEIGHT SET

Metric Weight Set, A nice set of quality weights are stored in a high impact open molded plastic hinged box. 1 gram and over are in bottle shaped, brass made and complete with forceps. All fractional weights are copper and aluminium made.



PA-92B

Weights set of 16 pcs, model PA-92B, each set includes 1 x 1gm, 2 x 2gm, 1 x 5gm, 2 x 10gm, 1 x 20gm, 1 x 50gm, Fractional weights include 1 x 10mg, 2 x 20mg, 1 x 50mg, 1 x 100mg, 2 x 200 mg, 1 x 500mg. Total 101.1 gram.



PA-92B-1

Weights set of 9 pcs, model PA-92B-1, each set includes 1 x 1gm, 2 x 2gm, 1 x 5gm, 2 x 10gm, 1 x 20gm, 1 x 50gm, 1 x 100gm. Total 200gram.



PA-92B-3

Weights set of 8 pcs, model PA-92B-3, each set includes 1 x 1gm, 2 x 2gm, 1 x 5gm, 2 x 10gm, 1 x 20gm, 1 x 50gm. Total 100gram.



PA-92B-4

Weights set of 10 pcs, model PA-92B-4, each set includes 1 x 1gm, 2 x 2gm, 1 x 5gm, 2 x 10gm, 1 x 20gm, 1 x 50gm, 1 x 100g and 1 x 200g. Total 400gram.



PA-92C-1

Weights set of 12 pcs, model PA-92C-1, each set includes 1 x 1gm, 2 x 2gm, 1 x 5gm, 2 x 10gm, 1 x 20gm, 1 x 50gm, 2 x 100gm, 1 x 200gm and 1 x 500gm. Total 1,000gm.



PA-92D-1

Weights set of 13 pcs, model PA-92D-1, each set includes 1 x 1gm, 2 x 2gm, 1 x 5gm, 2 x 10gm, 1 x 20gm, 1 x 50gm, 2 x 100gm, 1 x 200gm, 1 x 500gm and 1 x 1000gm. Total 2,000gm.



PA-92B-2

Economy Weight Set of 8 pcs, Model PA-92B-2, includes 1 x 1gm, 2 x 2gm, 1 x 5gm, 2 x 10gm, 1 x 20gm and 1 x 50gm weight. (without forceps), supplied in plastic bag. Total 100gm.