



PA-122C

“Globe” Electronic Stopwatch, Model PA-122C. A precise inexpensive quality instrument with accuracy to 0.01 second and with a 60 minute digital readout. 6 digit and 6 functions (hour, minute, second, month, date and day). Complete with battery, lanyard and instruction sheet with CE marked.



PA-122-6

Stop Clock, Model PA-122-6, A mechanical stop clock for all kinds of timing experiments, and is an educational instruments necessary for middle schools, with start, stop and reset buttons mounted on the top of the sturdy steel case, diameter of dial is 80mm divided 0-60 x 0.2 sec., and registers to 60 minutes. The clock will run for 24 hours on a single winding.



PA-122-1 & PA-122-2

Model PA-122-1
1/5 second and registers to 30 minutes.

Model PA-122-2
1/10 second and registers to 15 minutes.



PA-122-3 & PA-122-4

Model PA-122-3
1/5 second and registers to 30 minutes.

Model PA-122-4
1/10 second and registers to 15 minutes.

Single Action Stopwatch, an excellent stopwatch that is ideal for student use. Operation is simple with start, stop, and reset controlled by successive crown depressions. Attractive metal case.

Interruption Type Stopwatch, This jeweled stopwatch can be started, stopped, and restarted by depressing the crown. A second side button resets the hand to zero. The easy-to-read numerals are red on a white background. Case is nickel-chrome-plated for long life.



PA-122-5

Mechanical Timer, Model PA-122-5, measures one minute intervals up to one hour. After setting the required time interval, a signal indicates the end of the timed period. A bell sounds when the timer has completed the selected segment.