

GS Series High Precision Balance



GS223(W)~623(W)



GS822/1202



GS2202~12001
GS2202W/4202W

Model	Standard	GS223	GS323	GS423	GS623	GS822	GS1202	GS2202	GS3202	GS4202	GS6202	GS8201	GS12001
	With built-in calibration weight	GS223W	---	GS423W	GS623W	---	---	GS2202W	---	GS4202W	---	---	---
Capacity		220g	320g	420g	620g	820g	1200g	2200g	3200g	4200g	6200g	8200g	12000g
Readability		0.001g	0.001g	0.001g	0.001g	0.01g	0.01g	0.01g	0.01g	0.01g	0.01g	0.1g	0.1g
Pan Size		Ø118mm, with windshield			142 x 170mm		160 x 180mm						

Mono-Metal Tuning –fork Sensor (MMTS)

“SHINKO” GS Series High Precision Balance equipped with MMTS which provides quick response, long-term stability and high resolution

ISO Standard

“SHINKO” GS Series High Precision Balance can meet ISO9001, JQA-2834 TSUKUBA standard



CT6002 High capacity electronic carat balance



Capacity: 6000 ct
Resolution: 0.01 ct
1200 x 0.002g

Feature:

- With windshield
- Pan size: Ø118mm

12 Selectable units

Weighing units:

G, ct, oz, lb, ozt, dwt, GN, tl (Hong Kong), tl (Taiwan), tl (Singapore, Malaysia), momme, tola
GS series have 3 selectable modes (weight, percentage and unit count)

5 preset channels

“SHINKO” GS series have 5 preset channel for preset 5 common used units. Units can be changed by simply pressing the switch.

CS603 Light & Compact Carat Balance



Capacity: 600ct,
Resolution: 0.001ct

Feature:

- Small and compact, size: 235 x 182 x 168mm (include windshield)
- With windshield
- Save place, elegant design windshield. Roll-top design allows frontal access makes weighing easy
- Pan size: Ø80mm

CB603 Light & Compact Carat Balance



Capacity: 600ct
Resolution: 0.001ct

Features:

- Dimension: 251 x 207 x 260mm (include windshield)
- With windshield
- Pan size: Ø80mm

RS-232 as standard feature for **GS Series** and Carat Balance Series



香港科學儀器社有限公司
Che Scientific Co. (Hong Kong) Ltd.

Hong Kong Head office:

Unit D, 18/F, Roxy Industrial Centre, 58-66 Tai Lin Pai Road, Kwai Chung, Hong Kong

TEL: (852) 2481 1323 FAX: (852) 2418 1302

E-mail: sales@chescientific.com Website: www.chescientific.com



浙江科通儀器有限公司
Zhejiang Fortune Scientific Co. Ltd.

TEL: (86 571) 8729 7500 FAX: (86 571) 8729 7511

E-mail: sales@fortunescientific.com.cn

Website: www.fortunescientific.com.cn

Guangzhou Branch:

TEL: (86 20) 8130 7388 FAX: (86 20) 8108 4925

E-mail: guangzhou@fortunescientific.com.cn