




“TALBOYS”牌 加热磁力搅拌器
“TALBOYS” Hotplate Stirrers

美国制造
Made in USA

Digital Hotplate-Stirrer 加热磁力搅拌器

Features 特点:

- Low profile design
薄型设计
- LED display shows set-point temperature
LED 显示屏显示设定温度
- Built-in memory allows users to recall last set temperature
内存记忆可记下最后一次设定的温度
- Microprocessor control
微电脑控制
- Robust heater allows heating surface to reach set point quickly.
高功率加热器使受热面快速达到设定值
- Ceramic or aluminum top plate materials available
可选陶瓷或铝制平台材料
- Built-in support rod holder with locking knob
内置支持杆插孔及锁紧旋钮
- Housing is made of a heat and chemical resistant polymer that is cool to the touch
外壳使用耐热和耐化学性物料
- Hot symbol warning light is illuminated when heat is turned on and remains on until top plate cools down.
附高温警示灯, 亮起时表示平台是高温, 不可接触



Specification 规格

Top plate dimensions 平台尺寸		4 x 4" (10.2 x 10.2cm)	7 x 7" (17.8 x 17.8cm)	10 x 10" (25.4 x 25.4cm)
Temperature range 温度范围	Ceramic plate 陶瓷平台	Ambient(室温) +5° to 500°C		
	Aluminum plate 铝干台	Ambient(室温) +5° to 400°C		
Speed range 搅拌速度		60 to 1600rpm		
Max Stirring capacity 最大搅拌容量 (H2O)		600ml	2500ml	6000ml
Overall dimensions 外部尺寸 (LxWxH)		27.4 x 16.8 x 10.9cm	37.6 x 25 x 10.9cm	45.5 x 33 x 10.9cm

Ordering Information 订购资料

Order code 订货编号	Plate dimensions 平台尺寸	Plate material 平台材料	Electrical 电源
984TA4CHSEUA	4 x 4" (10.2 x 10.2cm)	Ceramic 陶瓷	230V 1.7 amps 400 watts
984TA4AHSEUA	4 x 4" (10.2 x 10.2cm)	Aluminum 铝	230V 1.7 amps 400 watts
984TA7CHSEUA	7 x 7" (17.8 x 17.8cm)	Ceramic 陶瓷	230V 4.6 amps 1050 watts
984TA7AHSEUA	7 x 7" (17.8 x 17.8cm)	Aluminum 铝	230V 4.6 amps 1050 watts
984TA0CHSEUA	10 x 10" (25.4 x 25.4cm)	Ceramic 陶瓷	230V 7 amps 1600 watts
984TA0AHSEUA	10 x 10" (25.4 x 25.4cm)	Aluminum 铝	230V 7 amps 1600 watts

Includes a power cable and a 3.8cm PTFE coated stir bar 连一条电源线和 3.8 厘米 PTFE 搅拌棒