

WEIFO®

“维和”牌 微电脑防潮箱
“WEIFO” Electronic Auto-Drying Cabinet

台湾制造
Made in Taiwan

Standard Series Electronic Auto-Drying Cabinet 实用型 微电脑防潮箱

数字显示湿度, 20-60%RH 湿度设定, DRY-205 型同时显示温度, 另有 Super Dry 型号, 10-60%RH 湿度设定



DRY-60

DRY-70

DRY-100

DRY-205

Economy Series Electronic Auto-Drying Cabinet 经济型 微电脑防潮箱

棒条式显示湿度, 30 / 40 / 50 / 60%RH 湿度设定



KD-70

KD-112

Model 型号	Standard Series 实用型				Economy Series 经济型	
	DRY-60	DRY-70	DRY-100	DRY-205	KD-70	KD-112
External Dim. (WxDxHmm) 外形尺寸 (阔×深×高 mm)	400x410x560	400x410x660	400x410x890	400x530x1280	400x410x560	400x410x890
Internal Dim. (WxDxHmm) 内部尺寸 (阔×深×高 mm)	398x355x480	398x355x580	398x355x820	398x475x1210	398x355x480	398x355x820
Capacity 有效容积	68L(公升)	82L(公升)	116L(公升)	228L(公升)	70L(公升)	116L(公升)
Door 门	SINGLE (单门)	SINGLE (单门)	DUAL(双门)	DUAL(双门)	SINGLE (单门)	SINGLE (单门)
Shelf 隔板	2 PCS (片)	2 PCS(片)	3 PCS(片)	5 PCS(片)	2 PCS(片)	3 PCS(片)
Weight 重量((约)	15kg (千克)	17kg (千克)	20kg (千克)	30kg (千克)	15kg (千克)	20kg (千克)

Professional Series Electronic Auto-Drying Cabinet 专业型 微电脑防潮箱

数字显示湿度及温 度, 20-60%RH 湿度设定, 另有 Super Dry 型号, 10-60%RH 湿度设定



DRY-420

DRY-610

DRY-530

DRY-1360

Professional Series 专业型						
Model 型号	DRY-300	DRY-420	DRY-520	DRY-610	DRY-530	DRY-1360
External Dim.(WxDxHmm) 外形尺寸(阔×深×高 mm)	890x480x930	1180x480x1020	1180x570x1020	1400x600x1055 连脚轮	920x660x1190 连脚轮	1200x730x1900 连脚轮
Internal Dim.(WxDxHmm) 内部尺寸(阔×深×高 mm)	888x425x850	1178x425x900	1178x515x900	1398x545x900	888x610x1070	1198x675x1740
Capacity 有效容积	320L(公升)	450L(公升)	545L(公升)	646L(公升)	578L(公升)	1402L(公升)
Door 门	DUAL (双门)	DUAL (双门)	DUAL (双门)	DUAL (双门)	DUAL (双门) / 另有 4 门可选 (前后双门)	6 PCS DOOR (六门)
Shelf 隔板	3 PCS (片)	3 PCS (片)	3 PCS (片)	3 PCS (片)	3 / 6 PCS (片)	5 PCS (片)
Weight 重量((约)	43kg (千克)	55kg (千克)	60kg (千克)	68kg (千克)	70kg (千克)	120kg (千克)

香港科学仪器社有限公司
Che Scientific Co. (Hong Kong) Ltd.
Tel: (852) 2481 1323
Fax: (852) 2418 1302

浙江科通仪器有限公司(绍兴)
Zhejiang Fortune Scientific Co. Ltd. (Shaoxing)
Tel: (86 575) 8812 0560
Fax: (86 575) 8812 0561

浙江科通仪器有限公司(广州)
Zhejiang Fortune Scientific Co. Ltd. (Guangzhou)
Tel: (86 20) 8130 7388
Fax: (86 20) 8108 4925

上海香科仪器贸易有限公司
Shanghai Xiang Ke Scientific Trading Co. Ltd.
Tel: (86 21) 6248 6170
Fax: (86 21) 6248 0647